



Check & Trace Salmonella 2.0

Save valuable lab time with our new molecular *Salmonella* serotyping workflow using real-time PCR multiplexing technology.

Get results in under two hours.



Get in touch with
our experts
www.checkandtrace.com



Rapid and actionable

Confirm *Salmonella* presence and identify the top 59 serotypes of concern in under 2 hours.

Typing is essential to eliminate *Salmonella* from the production chain.

Multiplex real-time PCR-based technology detects 27 target genes in parallel.



Efficient and easy-to-use workflow

Easy sample preparation process.

Plug-and-play assay with ready-to-use reagents saving precious lab time.

Use standard laboratory equipment and a Bio-Rad real-time PCR cycler to run the test.

Access and manage data anytime, anywhere with our cloud-enabled services.



Objective results

Proprietary software interprets results and provides a simple readout of the *Salmonella* serotype, based on a unique profile of genetic markers.

How it works: Check & Trace Salmonella 2.0 workflow

1. Sample preparation	2. Amplification	3. Result interpretation
<ul style="list-style-type: none"> – Take sample from pure culture (XLD or nutrient agar medium plate) using colony sampler – Add buffer, mix and heat to extract DNA 	<ul style="list-style-type: none"> – Mix DNA sample with Master Mix (reaction mix) – Aliquot reaction mix into prefilled strips – Run real-time PCR (90 mins)* 	<ul style="list-style-type: none"> – Upload cycler results* to the Check & Trace analysis software to receive results in minutes

*For use with Bio-Rad CFX96 / CFX Opus 96

Salmonella serotypes covered

Salmonella discrimination is of paramount importance when carrying out source-tracking investigations. Early *Salmonella* detection and discrimination with Check & Trace Salmonella 2.0 (Cat. No. 18-2020) enables rapid decisions and responses.

Abaetetuba	Choleraesuis	Havana	Molade	Paratyphi B. (Java)	Tennessee
Agona	Corvallis	Heidelberg	Montevideo	Poona	Thompson
Alachua	Cubana	Idikan	Muenchen	Pullorum	Typhimurium
Albany	Derby	Infantis	Muenster	Reading	Typhimurium monophasic (1,4,[5],12:i:-)
Bovismorbificans	Dublin	Javiana	Newport	Rissen	Uganda
Braenderup	Enteritidis	Kentucky	Ohio	Saintpaul	Virchow
Brandenburg	Gallinarum	Livingstone	Oranienburg	Sandiego	Worthington
Braenderup	Give	London	Orion	Schwarzengrund	Yuruba
Bredeney	Goldcoast	Mbandaka	Ouakam	Senftenberg	4,[5],12:d:-
Cerro	Hadar	Minnesota	Panama	Stanley	

Contact information

Visit: www.checkandtrace.com
 Support & orders: info@check-points.com
 Phone: +31 (0) 317 453 908



INDICAL
 BIOSCIENCE

Check-Points B.V.
 An INDICAL BIOSCIENCE company
 Binnenhaven 5, 6709 PD Wageningen,
 The Netherlands

Please carefully read eIFU for intended use, performance characteristics and product limitations. Trademarks: Bio-Rad®, CFX96™ (Bio-Rad Laboratories, Inc.); INDICAL® (INDICAL BIOSCIENCE GmbH). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-CTS2-EN-230113
 © 2023 Check-Points B.V.